
**Varian Gemini2000 (300 MHz) NMR Check-Out Procedure
for Year 2000
May 15th—June 3rd, 2000**

1. Become familiar with the operation of the instrument with an authorized user.
(An authorized user must be present at all times while you are training)
2. Contact Tanase (tel 20-3399, tanase@cc.nara-wu.ac.jp) to schedule a check-out.
3. Demonstrate at the check-out that you can run a 1D ¹H spectrum of a test sample and that you have a basic familiarity with the following procedures.
You can see the manual by Dr. Mikata and related ones.

Check-Out List

- 1) Inserting and spinning the sample.
- 2) Loading the correct parameter file for your sample.
- 3) Locking and shimming:
 - a) Centering the lock signal and locking
 - b) Adjusting lock phase
 - c) Optimizing the Z1 and Z2 zhims.
 - d) Loading a standard shim file.
 - etc
- 4) Modifying the acquisition parameters:
 - a) be familiar with pulse sequence
 - b) change acquisition time (why and how)
 - c) change pulse width (why and how)
 - d) change the relaxation delay (why and how)
 - e) change observation window (why and how)
 - etc
- 5) Obtaining FID and spectrum
 - a) adjusting the receiver gain (why and how)
 - b) check FID data (why and how)
 - c) weighting Fourier transformation
 - d) phasing the spectrum
 - e) setting the chemical shift reference
 - f) integrating the peaks
 - g) plotting
 - etc
- 6) Saving and loading the FID data.
- 7) Ejecting the sample and setting up the machine for a next user.

Good luck!!
Tomoaki Tanase (Room C323)
tel 0742-20-3399
tanase@cc.nara-wu.ac.jp

What you can do as routine measurements with Gemini2000

1D	1H	Nondecoupling (standard)	o <recommended>
		Homo Spin Decoupling	o
		Homo Gated Doupling	o
	13C	Broad Band Decoupling (standard)	o <recommended>
		Offresonance Decoupling	o (not good)
		Gated Decoupling	o
		Inverse Gated Decoupling	o
		Attached Proton Test (like Dept135)	o <recommended>
	31P	DEPT	o(not good)
2D	31P	Inverse Gated Decoupling (standard)	o
		19F	Nondecoupling (standard)
			not checked
	cosy	HH Cosy(AV)	o <recommended>
	dqcosy	Double Quantum Trans. HH Cosy (PS)	o
	noesy	HH Noesy(PS)	o
	cosyps	HH Cosy(PS)	not checked
	hetcor	CH Cosy(AV)	o (not good)
VT	hmqc	HC	not checked
	hmbc	HC long range correlation	x
	high temp	vttype=0 5~100 °C	o
	low temp	vttype=2 -80 ~ 5 °C	o

2000.5 by TT