

Instruction Manual to input reagents data into the database file "Reagent List" prepared by File Maker Pro 4.0 (for Macintosh)

All reagents in TT lab must be filed !!

The reagent database can be used as Material Safety Data

Please click the note yellow marks.

Do NOT delete this card even when the bottle is empty. In this case
please input a comment in the Notes field.

REAGENT		Nara Women's University, Department of Chemistry Tanase Group	
Name	chloromethyltrimethylsilane, ClCH ₂ Si(CH ₃) ₃ 		
Formula	C ₄ H ₁₁ ClSi 		
Notes	<div>毒劇物 クロロメチルトリメチルシラン 125 g * 1 冷暗所保存、火気厳禁、毒劇物 dpmp合成用 </div> <div>(注) 化合物名は正式名称以外にも慣用名、示性式等を、で区切って入力する。行の先頭に数字やcis,trans等の記号が来ないようにする。組成式はCHNOアルファベットの順で元素を書き1個の場合にも1を記入する。例えば、C₆H₁₀N₂O₃Co1Zn1。</div> <div>Class ORG, Si Place FRZ2  Date 99/01/05</div>		
Caution!	<div><ul style="list-style-type: none">● 発火性あり● 毒性あり● 引火性あり● 悪臭あり● 空気に不安定● 湿気に不安定● 光に不安定● 熱に不安定● 安定</div> <div>Handling! <ul style="list-style-type: none">● 厳重窒素下● 低温下● 厳重遮光● デシケーター● ボンベ● その他...</div>		